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BioPAT® Xgas

English - page 3

In cases involving questions of interpretation, the German-language version shall prevail.

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The enclosed CD contains the operating instructions as a PDF file in various international languages.

If the CD is missing, you can obtain a copy form us by specifying the order number:

Order number: 85037-548-53

Sartorius Stedim Biotech GmbH Technical Editorial Department August-Spindler-Strasse 11 37079 Goettingen, Germany

tech.pubs@sartorius-stedim.com www.sartorius-stedim.com

System Requirements:

- Windows, MacOS X
- Browser with JavaScript enabled
- PDF-Reader

Use start.html to run the application

Important Note:

The data and information contained in this manual were compiled with the greatest of care. Notwithstanding the scrupulous attention paid to preparing this document, no absolute guarantee of its accuracy can be given. Please notify us if any information is missing that is important for you or if you find any errors.

Contents

1.	Intro	oduction 4
	1.1	Function 4
	1.2	Target Group 4
	1.3	Notes on Using these
		Instructions 4
2.	Safe	ety Instructions 6
	2.1	General Information 6
	2.2	Authorized Personnel 6
	2.3	Intended Usage 6
	2.4	
	2.5	Process-related Hazards 7
3.	Desi	gn and Function 8
	3.1	
		3.1.1 Pressure Measurement . 10
		3.1.2 Humidity Measurement 10
		3.1.3 CO ₂ Measurement 10
		3.1.4 O ₂ Measurement 10
4.	Insta	allation11
	4.1	Equipment11
		4.1.1 Unboxing and
		Assembling 12
	4.2	Setup
	4.3	Assembly
	4.4	Connection
	4.5	Arrangement on the
		Bioreactor and Assembly 22

5.	Ope	rating	Information	23
	5.1	Safet	y Instructions	23
	5.2	In-pro	ocess Exhaust	
		Air Ar	nalysis	23
		5.2.1	Switching on the	
			BioPAT® Xgas	24
		5.2.2	In-process	
			Measurement	25
		5.2.3	Sensor Status	
			Display via LEDs	
	5.3		leshooting	28
		5.3.1	Faults or Malfunctions	
			from Moisture	
			in the Exhaust Air	
		5.3.2	Leaktightness Test	29
6.	Maiı		nce and Service	30
	6.1		tenance	
			e Customer	
			Calibration	
		6.1.2	Filter Replacement	31
7.	Арр	endix.		34
	7.1	Calib	ration Table	34
8.	Tech	nical	Data	35

1. Introduction

1.1 Function

This operating manual describes the safe installation and use of the BioPAT® Xgas. The operating personnel must read these instructions prior to operational startup. These instructions should always be available at the workplace.

1.2 Target Group

These operating instructions are designed for qualified personnel. Only qualified personnel may set up BioPAT® Xgas, connect it to the bioreactor or perform in-process operation or maintenance.

1.3 Notes on Using these Instructions

All information in these operating instructions considers the state of the art, the applicable standards and guidelines, and our many years of experience and knowledge.

These operating instructions provide you with all the information necessary to install and operate the BioPAT® Xgas (referred to as "device" in the following).

The operating instructions must be read and used by all persons entrusted with the operation, service, cleaning, and troubleshooting of the device. This applies particularly to the safety instructions listed.

- ▶ Please read these operating instructions carefully and completely before using this device.
- ► These instructions are part of the product. Keep them in a safe and easily accessible place near the device's site of installation.
- ► If the instructions should be lost or misplaced, please contact Sartorius for a replacement or download the latest version from our website: www.sartorius.com

Symbols Used

As a means of instruction and direct warning against hazards, all text statements to be particularly noted in these operating instructions will be marked as follows:

This symbol denotes a possible danger with risk that moderate or minor injury may result if it is not avoided

This symbol denotes a danger with low risk that could result in property damage if it is not avoided.

This symbol indicates useful information or a function setting on the device, or that caution should be exercised while working.

The following presentations will also be used:

- Texts that follow this mark are lists.
- ► Texts that follow this mark describe activities that must be carried out in the specified order.
- Texts that follow this mark describe the result of an action.
- " Texts in quotes are references to other chapters or sections.

2. Safety Instructions

2.1 General Information

Sartorius Stedim Biotech GmbH supplies the BioPAT® Xgas after careful testing in ready-for-operation condition.

This manual contains instructions for the safe operation of the BioPAT® Xgas. You may only operate the device under conditions that comply with the specifications stated in this manual and the information on the manufacturer's ID label.

2.2 Authorized Personnel

The following persons may carry out the actions described here:

Operations	Persons
Assembly and connection at to the bioreactor; function test, basic calibration during servicing and maintenance	Trained and qualified personnel of the operator or service technicians of Sartorius Stedim Biotech GmbH
Operation, in-process 1-point calibration	Personnel who has read this documentation and is familiar with the operation of bioreactors
Factory maintenance and service	Service Department of Sartorius Stedim Biotech GmbH or the sensor manufacturer.

2.3 Intended Usage

The BioPAT® Xgas is a gas sensor for measuring carbon dioxide and oxygen contents in the defined concentration range and under the conditions described in the

▶ Technical Data. It is used to monitor metabolic activities during biological processes such as fermentations.

2.4 Limitations of Use and Misuse

Danger of malfunctions and damage to the sensor under impermissible measuring conditions. The BioPAT° Xgas is not designed according to the ATEX Directive.

Do not use the sensor

- in flammable or explosive atmospheres or with such gases!
- in gases with polymers or silicone components
- in gases with halogens like F, Cl, Br etc., CFC
- in gases with SOx and H₂S.

Prevent water or humidity (RH <75%) from entering into the sensor when the sensor is switched on. Dry the sensor when humidity condenses in the exhaust line and condensate has entered into the sensor.

The BioPAT® Xgas is unsuitable for use as safety equipment for system gas monitoring or as a gas warning device. This constitutes misuse.

2.5 Process-related Hazards

Risk of poisoning for persons when CO_2 is released into the room. Risk of inflammation for flammable materials if gases enriched with oxygen escape into the room air.

Determine the max. concentrations of gas released into the room and evaluate whether they could lead to hazardous situations, e.g. if the room ventilation fails.

Ensure that your work area is well ventilated or connect the BioPAT* Xgas to the available exhaust air fittings in your laboratory.

Check all connections for leaks.

3. Design and Function

The BioPAT® Xgas combines one oxygen and one carbon dioxide sensor featuring automatic compensation of humidity and pressure into one compact analyzer. The gas to be analyzed flows through the integrated flow adapter past the measuring chambers. During this process, the pressure, humidity, oxygen and carbon dioxide concentration are measured.

The sensors measure the gases in the concentration ranges and under the conditions for which they were designed; observe the ➤ Technical Data. The analyzer transfers the measurements to the control unit of the bioreactor, e.g. the DCU systems of bioreactors BIOSTAT® Cplus, B and B-DCU II.

The measurement and control system shows measurements on the display, as configured in the settings. Observe the operating instructions for the measurement and control system.

3.1 Measurement Principles

The analyzer features 4 different measuring components:

- 1. Pressure measurement
- 2. Humidity measurement
- 3. CO₂ measurement
- 4. 0, measurement


Fig. 1: schematic illustration of complete installation

In order to prevent humidity contained in the exhaust air in the BioPAT® Xgas from condensing, the sensor housing is heated. Warming up prior to initial startup as well as every time the device is disconnected from the power supply will take about an hour. The sensor won't provide any valid readings during heating.

3.1.1 Pressure Measurement

A piezoresistive silicon pressure sensor is used to measure the pressure.

3.1.2 Humidity Measurement

A capacitive polymer sensor is used to measure the humidity.

3.1.3 CO₂ Measurement

The CO₂ sensor comprises an IR radiation source, a detector and the reflection measuring cell (Fig. 2).

The infrared beam is reflected by the gas-filled measuring adapter. The light, now weakened by the CO₂, is measured by the detector.

A translucent sapphire disc prevents the optical components from getting soiled by constituents of exhaust air.


Fig. 2: Measuring setup of the CO₂ sensor

Fig. 3: Measuring setup of the O_2 sensor

3.1.4 O_2 Measurement

The BioPAT® Xgas oxygen sensor is based on an oxygen pump cell and runs independently of reference gases. Applying voltage to the cell moves the oxygen ions from the cathode to the anode and generates current. Through a diffusion barrier in front of the cathode, a saturation current is generated by increasing the supply voltage, which allows measurement of the ambient oxygen concentration.

4. Installation

4.1 Equipment

Danger of malfunctions when damaged during transport and when operated under process conditions that do not comply with the device specifications.

When delivered for or with bench-top bioreactor systems, e.g. BIOSTAT® B or B-DCU II, the BioPAT® Xgas and accessories are packaged in a safe-to-transport manner.

When delivered with sterilizable-in-place (SIP) bioreactors, e.g., BIOSTAT® CPlus, the BioPAT® Xgas sensor is installed on the bioreactor frame. The sensor accessories are packaged separately in a safe-to-transport manner.

- 1. Carefully unpack the packaged devices and the connecting accessories.
- 2. Check if the delivery is complete.
- Check all delivered components for any damage. If transport damage has occurred, send a claim as soon as possible.

- The analyzer BioPAT® Xgas.
- Flow adapters for connection to the exhaust line.
- For bench-top bioreactors: connection hoses, e.g., silicon hoses with diameter
 6 mm, for connecting the exhaust filter to the BioPAT® Xgas Sensor.
- For sterilizable-in-place (SIP) bioreactors: installation kit already attached to the bioreactor frame.
- Signal transmission cable to the BioPAT® Xgas on the control unit.

Fig. 4: Analyzer, signal cable and accessories for mains connection (A) Spare kit 1: SB-20-02-0002, Power supply for BioPAT® Xgas

- When delivered with the external power supply, the male connector must match the power connection in laboratory. The delivery contains the matching plug adapter.
- A replacement BioPAT® Xgas power supply kit pack can be ordered from Sartorius (SB-20-02-0002), see "Fig. 5: Cable set included".

Fig. 4a: Unboxing


Fig. 4b: Preparing the assembling (B) Spare kit 2: SB-20-02-0001, Spare part kit for BioPAT® Xgas


Fig. 4c: Filters

4.1.1 Unboxing and assembling

The flow adapter of the BioPAT° Xgas is not preassembled and will be delivered in a small box together with the BioPAT° Xgas (fig 4a). Before the installation, it must be connected to the BioPAT° Xgas.

For the assembling you will need (fig. 4b):

- 1. The BioPAT® Xgas (A)
- 2. The TORX T20 screwdriver (D)
- 3. The mounting base plate (C)

the filter set (fig. 4c).

- 4. The flow adapter (E)
- 5. Bag with 2 flat head screws M4x6 TORX for the wall mounting and 4 screws M4x16 TORX for the flow adapter (B).
- 6. Bag with flow path filter set BioPAT® Xgas (B).

To assemble the BioPAT® Xgas proceed in the following steps:

. First, place the filters and gaskets on the flow adapter. All required filters and gaskets are in the filter set for the BioPAT® Xgas (SB-20-02-0001, Spare part kit for BioPAT® Xgas). For each recess there are two different filters.

There is a flat filter (fig. 4c: A) and a thicker filter (fig. 4c: B). The flat filter is hydrophobic and the thicker is the particulate filter. For each recess

there is a gasket, a flat filter and a thicker filter in

Fig. 4d: Filters and gaskets

Fig. 4e: Insert flat filter

Fig. 4f: Place gaskets

Fig. 4g: Insert second filter

Fig. 4h: Fit in filter

- 2. Insert the flat filter into the recesses (fig. 4e).
- 3. Place the gaskets in the recesses (fig. 4f)
- 4. Insert the thicker filters into the recesses (fig. 4g). If the filters don't fit in the recesses, it might be necessary to push them gently into place (fig. 4h).

Fig. 4i: Put screws in

Fig. 4j: Tighten screws

Fig. 4k: Fit base plate

- Now place the 4 TORX screws M4x16 (fig. 4b: E) in the bores of the BioPAT® Xgas without screwing. This will make it easier to place the BioPAT® Xgason the right spot on the flow adapter (fig. 4i).
- 6. Place the BioPAT® Xgas on the flow adapter. Tighten the 4 screws with the screw driver TORX.
- 7. Put the wall mounting plate on the flow adapter and use the two flat head screws M4x6 (fig. 4b: E) to fix the base plate to the sensor (fig. 4k).

4.2 Setup

Danger of damage to the control and supply unit of the bioreactor.

The BioPAT° Xgas module weighs approx. 3 kg. Do not place the module on the control and supply unit of the bioreactor and do not fasten it to its side panels.

The housing panels cannot carry the weight of the BioPAT® Xqas.

4.3 Assembly

When using bioreactor bench-top systems, mount the BioPAT° Xgas module on a flat surface with non-slip rubber feet.

- Observe the lengths of the cables and connection lines.
- 2. Place the BioPAT® Xgas on the work table next to the bioreactor so that you can comfortably connect all equipment.

Mounting the BioPAT Xgas on a fixed vertical surface is possible with the approriate 6.6 mm screws.

When using sterilizable-in-place (SIP) bioreactors, the BioPAT® Xgas module is attached to the bioreactor frame using an installation kit.

- Avoid contact between the signal cable and the outer surface of the vessel.
- 2. Ensure that no hot steam is directed at the sensor and signal cable.

4.4 Connection

Making connections while the BioPAT® Xgas is still connected to power or with the wrong kind of voltage can cause damage.

Use only cables from the equipment supplied with the BioPAT® Xgas and | or the bioreactor equipment. The module is not equipped with an ON | OFF switch. First, leave the control unit and | or supply unit of the bioreactors switched off and | or disconnect the external power adapter from the mains current. First connect the cables for the signal transmission and power supply.

Connect the exhaust line from the BIOSTAT® B and B-DCU II

- Make sure that the lines of the culture vessel contain no condensation.
- 2. Connect the exhaust line to the BioPAT® Xgas.
- 3. If exhaust air entering the room air can be regarded as safe, you can leave the BioPAT® Xgas outlet open.

If the exhaust air contains gases in potentially hazardous concentrations, connect the outlet to the exhaust air fittings in the laboratory.

Fig. 5: Cable set included

Connecting the Cables to the BIOSTAT® B and B-DCU II

 Use only the cables that are especially designed for this purpose. Observe the specifications below and the "Technical Documentation on the Bioreactor":

Wait until the BioPAT® Xgas has reached room temperature.

Connecting the Y cable

Take the Y cable of the BioPAT® Xgas.

Connect the female connector with the three cable outputs to the left plug of the BioPAT® Xgas.

Connect the BioPAT® Xgas to the power supply.

Connecting to BIOSTAT® B (Serial No. ≤ 1000)

Connect the analog signal outputs of the BioPAT° Xgas to the analog inputs on the Control Tower of the BIOSTAT° B as follows:

Take the analog Y signal cable.

Connect the plug of the Y cable labeled "Xgas" to the right female connector of the BioPAT® Xgas.

Connect the plug of the Y cable labeled "Ext. Sig.-A" to the female connector "Ext.Sig.-A1/2" of the BIOSTAT® B.

Connect the plug of the Y cable labeled "Ext. Sig.-A" to the female connector "Ext.Sig.-B1/2" of the BIOSTAT® B.

Connecting to BIOSTAT® B-DCU II

Connect the analog signal outputs of the BioPAT° Xgas to the analog inputs on the Control Tower of the BIOSTAT° B-DCU II as follows:

Take the analog Y signal cable.

Connect the plug of the Y cable labeled "Xgas" to the right female connector of the BioPAT® Xgas.

Connect the plug of the Y cable labeled "Ext. Sig.-A" to the female connector "Ext.Sig.-A" of the BIOSTAT® B-DCU II.

Connect the plug labeled "Ext. Sig.-B" of the Y cable to the female connector "Ext.Sig.-B" of the BIOSTAT® B-DCU II.

Connecting to BIOSTAT® B (Serial No. > 1000)

Connect the analog signal outputs of the BioPAT® Xgas to the analog inputs on the control tower of the BIOSTAT® B as follows:

Take the analog signal cable.

Connect the cable plug labeled "Xgas" to the right-hand socket of the BioPAT® Xgas.

Connect the cable plug labeled "Ext. Sig.A1/B1" or "Ext. Sig.C1/D1" to the "Ext.Sig.-A1/B1" or "Ext.Sig.-C1/D1" socket of the BIOSTAT® B depending on the DCU signal configuration for off-gas analysis.

Connecting to BIOSTAT® Cplus

Connect the analog signal outputs of the BioPAT® Xgas to the control tower of the BIOSTAT® Cplus as follows:

Take the analog serial connection cable.


Connect the cable plug labeled "Xgas" to the right-hand socket of the BioPAT® Xgas.

Connect the cable plug labeled "Ext. Sig.-A" to the "Ext. Sig.1-2" socket of the BIOSTAT® Cplus.

4.5 Arrangement on the Bioreactor and Assembly

Fitting the BioPAT® Xgas to a BIOSTAT® Cplus:

A: BioPAT® Xgas Sensor & flowpath module

B: Support

C1/C2: Valves for switching the exhaust feed
D: Pneumatic pressure control valve
E: Pipe to BioPAT® Xgas sensor

F: Pipe between exhaust feed and pressure control valve. G1/G2: Outlets into the ambient air or for connection to the

laboratory vent, if necessary.

5. Operating Information

5.1 Safety Instructions

Hazards for users from the composition of the exhaust air:

- CO₂: Danger of suffocation from release into the room air:
- O₂: Combustion of flammable materials at high
 O₂ concentrations or potential production of oxy-hydrogen gas with H₂.

When channeling into room air, ensure that your work area is well ventilated. If constituents, e.g. generated by reactions with other substances, may be hazardous, channel the exhaust air into ventilation fittings on the laboratory side. Temperatures in the exhaust air > 40°C or oxygen-free exhaust air can damage the sensors in the BioPAT° Xgas. The pure nitrogen used during zero-point calibration of the DO sensor in the culture vessel must not escape into the BioPAT° Xqas.

At process temperatures >40°C and during zero point calibration of the DO sensor with nitrogen, leave the BioPAT° Xgas disconnected from the exhaust line and | or block the inlet lines or disconnect them temporarily.

5.2 In-process Exhaust Air Analysis

Danger of burns! After connecting the power supply, the housing of the BioPAT® Xgas heats up due to warming. Do not touch it any longer or only when wearing gloves.

Compare the process conditions with the bioreactor specifications in the "Technical Data". Do not connect the device for gases or under process and ambient conditions it is not designed for.

5.2.1 Switching on the BioPAT® Xgas

After initial start-up and when the power supply is disconnected, the BioPAT° Xgas requires a warming phase of about 1 h until it can operate within the limits of the measuring range, see "Technical Data".

LEDs on the device display the operation status. Displayed measurements indicate the zero current conducted during warm-up, not the actually measured gas concentrations.

- 1. The BioPAT® Xgas is activated when the bioreactor is switched on or when the power supply is connected. Wait for the device to warm up (approx. 1 h). The LED rings around the buttons light up and indicate the sensor status.

 During the warm-up phase, the sensor supplies a current signal of approx. 2.3 mA
- 2. Perform 1-point calibration during initial startup and always when the analyzer indicates you do so. Disconnect the connection or lines from the bioreactor (exhaust air filter) and channel fresh air (0.4% CO₂, 20.97% O₂ by volume, residual N₂) at a flow rate of min. 200 ml/min. Observe the specifications in the "Data Sheet". Flushing air with a different composition will falsify later measurements.
- Press the two buttons on the BioPAT® Xgas for 5 seconds.

As soon as the left LED is lit up continuously in green, the analyzer is ready for operation.

Fig. 6: 1-point calibration

For autoclavable bioreactors:

When disconnecting from the exhaust line during preparation of the bioreactor for the process, e.g. for autoclaving, you must not switch off the BioPAT® Xgas. The displayed measurement signal then indicates the measurement of the ambient air

- On bioreactors with autoclavable culture vessels, remove any exhaust line already connected to the BioPAT® Xgas.
- 2. Prepare the bioreactor for the process [➤ Operating manual for the bioreactor], e.g. sterilization, calibration of the DO sensor and inoculation.

5.2.2 In-process Measurement

Temperatures in the exhaust air >40°C or oxygen-free exhaust air can damage the sensors in the BioPAT° Xgas. The pure nitrogen used during zero-point calibration of the DO sensor in the culture vessel must not escape into the BioPAT° Xgas. Humidity and condensation in the exhaust air can damage the sensors.

Ensure that exhaust cooling is always active during measurement.

Monitor whether condensation forms in the exhaust line. Abort the measurement if condensation enters the BioPAT® Xgas. Dry the BioPAT® Xgas as described below [> Troubleshooting].

Bioreactors BIOSTAT® B or B-DCU II

- 1. After autoclaving, set up the culture vessels at the work place. If required, calibrate the zero point of the built-in DO sensor with pure nitrogen.
- 2. Ensure that the process temperature is < 40 °C.
- 3. Connect the exhaust line to the BioPAT® Xgas.

BIOSTAT® Cplus bioreactors

- If the zero point of the DO sensor with oxygen is calibrated, switch the exhaust line to the bypass via the solenoid valves.
- 2. Calibrate the zero point of the built-in DO sensor with nitrogen.
- 3. Ensure that the process temperature is below 40°C.
- 4. Switch the exhaust line to the BioPAT® Xgas via the solenoid valves. For bioreactors with automatic valves, check that the valves have switched using the control.

If the exhaust is not directed to the ambient air, ensure that both the BioPAT® Xgas and the bypass are connected to their joint line.

Conducting measurements with the BioPAT® Xgas

- Every time the device is connected to or disconnected from the power supply, it will take about an hour to warm up. The displayed measurement signal is related to the current signal of approx.
 2.3 mA from the sensor during the warm-up phase
- If the measurement and control system of the bioreactor only indicates the measured values, monitor the course of process on the operator display.
- If the measurement and control system of the bioreactor can log the measured values, configure this log as described in the [▷ Operating manual for the control system].

5.2.3 Sensor Status Display via LEDs

LED status	Description	Possible Causes
green	Sensor is ready for operation	
yellow	Warm-up phase, can take up to 60 minutes	The sensors must reach operating temperature, otherwise the measurements are incorrect
green yellow	Prompt to perform a 1-point calibration	Always at the initial startup, or when the sensor signal shows extreme changes
green flashing	1-point calibration running	Calibration started by pressing the two buttons for 5 seconds
green yellow flashing	Factory calibration required	Strong signal alteration due to age of the sensor or incorrect 1-point calibration
red	Sensor not working	Sensor currently starting up; signal too weak; sensor defective


Fig. 7: Buttons with LEDs

5.3 Troubleshooting

5.3.1 Faults or Malfunctions from Moisture in the Exhaust Air

Humidity, e.g. condensation in the exhaust line, can damage the sensors in the BioPAT® Xgas or disrupt the exhaust measurement.

In the event of moist exhaust air that is warmer than the ambient air, condensation can form in the exhaust line and be forced into the module. Check the exhaust line regularly for any condensation that may have formed.

Allow the condensation to drain off and dry the BioPAT® Xgas with compressed air.

You do not need to interrupt the measurement. In doing so, note that the device is now determining the gas concentrations in the air being pumped through.

System interruptions on BIOSTAT® B or B-DCU II

- 1. Disconnect the tubing from the inlet to the BioPAT® Xgas.
 - Allow the condensation to drain off, e.g. into a vessel or down the lab drain.
- 2. Connect clean, dry compressed air and flush the BioPAT® Xgas.
- 3. Reconnect the exhaust air segment tubing to the BioPAT® Xgas.

If larger amounts of condensation occur, and you continually have to dry the BioPAT® Xgas, you can channel the exhaust air through a condensation trap [> Information on request].

5.3.2 Leaktightness Test

Ambient air entering the exhaust air disrupts the measurement. Check whether the tubing, lines, connections and the BioPAT® Xgas are leaktight

- 1. Block the gas outlet of the BioPAT® Xgas or clamp off the lines to the exhaust air fittings in the laboratory.
- 2. Check whether gas is escaping from tubing, lines, connections or the housing.
- Replace defective tubing, tighten loose screws and | or fasten the tubing with hose clamps.
 On pipelines with Tri-Clamp connections, check and replace the Tri-Clamp seals and gaskets.
- If gas is escaping out of the housing of the BioPAT® Xgas, check that the flow adapters are assembled leaktight, e.g. when you have maintained the filter installed there [> Maintenance Instructions]

6. Maintenance and Service

6.1 Maintenance by the Customer You must calibrate the BioPAT® Xgas sensors regularly. Observe the LED displays on the BioPAT® Xgas analyzer:

6.1.1 Calibration

- Perform the 1-point calibration during initial startup and always when the analyzer indicates you do so.
- Aerate once monthly for approx. 30 minutes with normal fresh air (20.97% by volume O₂ | 0.04% by volume CO₂) as per specification (see Data Sheet). Next, press the two keys at the same time for 5 seconds (Fig. 8).
- For the recommended annual recalibration, return the sensor to Sartorius Stedim Systems GmbH.


Fig. 8: Screws for attaching the sensors to the flow adapter

6.1.2 Filter Replacement

The flow adapter contains a coarse filter, which you have to maintain and replenish whenever there are impurities or obstructions of the gas flow. Spare parts are obtained on request. For filter replacements:

- Disconnect all connected lines from the BioPAT°
 Xgas or disassemble the module from the bioreactor. Remove the mounting plate as necessary.
- 2. Unscrew the four screws visible from top (see Fig. 8). Remove the base plate and place it to the side so that the filters are visible (Fig. 8)
- 3. Remove the filters from the recesses and insert new filters.
- 4. Replace O-rings with pressure marks, porous surfaces or hairline cracks.
- Re-fasten the base plate and install the BioPAT® Xgas on the bioreactor and | or mount and connect.

Replace the BioPAT® Xgas if the sensors are defective.

An LED lighting up continuously in red after the startup and warm-up phase indicates malfunctioning or defective sensors.

On-site repairs are not possible. Inform the Service Department of Sartorius Stedim Biotech GmbH.

Sartorius Stedim Biotech GmbH recommends having the sensor checked and recalibrated once a year at the factory. Moreover, send the sensor back in the event of an incorrect or failed 1-point calibration, in the event of strong signal alteration without major changes to the measuring conditions.

Sartorius Stedim Biotech GmbH shall not be liable for repairs made by the customer or any resulting damage, especially in the event of

- use of non-approved spare parts.
- non-approved technical changes.
- Inspections can be carried out on-site by authorized service personnel or by any service representative of Sartorius Stedim Biotech GmbH.
- 2. In the event of service requests or warranty claims, please notify your representative of Sartorius Stedim Biotech GmbH or contact:

Sartorius Stedim Systems GmbH Service Robert-Bosch-Strasse 5–7 D-34302 Guxhagen, Federal Republic of Germany Phone +49 (0) 5665-407-0 / Fax +49 (0) 5665-407-2200

3. You can return defective devices or parts.

Returned devices must be clean and in hygienically flawless condition and packed carefully, e.g. in accordance with the safety guidelines applicable to the process.

- Inform the Service Department in the event of a return.
- 2. Provide a Declaration of Decontamination with the return shipment.

The Service Department must be able to read it before opening the packaging. The information must detail the media and microorganisms the **BioPAT**[®] **Xgas** has come into contact with and the measures taken to disinfect and decontaminate it.

7. Appendix

7.1 Calibration Table

Only have the complete calibration conducted by Sartorius Stedim Biotech GmbH Service.

Conduct the monthly 1-point calibration as described in Section ▶ "6.1 1-Point Calibration" on page 27.

Also calibrate this way every time after separating sensor head and measuring adapter.

Log the calibrations performed, e.g. in the table below.

Date	Procedure	Conditions	Name	Signature
	Complete calibration	25 °C, 1 bar	BioPAT® Xgas	
	1-Point calibration			

8. Technical Data

O₂ Sensor Unit

Gas	0_2
Concentration range	1 – 50% by volume, others available upon request
Measuring principle	Zirconium dioxide
Accuracy	< 0.2% FS2 ± 3% Rdg2 on value
Drift	< ± 2% value year
Lifetime of sensor element	Approx. 15,000 operating hours approx. 2 years

Temperature inside the sensor unit 580°C | 1076°F

CO₂ Sensor Unit

Gas	CO ₂
Concentration range	0 – 10% by volume, others available upon request
Measuring principle	Infrared
Accuracy	< 0.2% FS2 ± 3% Rdg2 on value
Drift	< ± 2% value year
Lifetime of sensor element	Approx. 25,000 operating hours approx. 3 years

Temperature inside the sensor unit 3°C | 5.4°F higher than process temperature

General Information

Temperature range	15 – 40°C 59°F – 104°F
Pressure range	0.8 – 1.3 bar 11.6 – 18.85 psi absolute pressure
Humidity range	0 – 100% relative operating humidity, integrated humidity compensation
Housing	Stainless steel, IP65
Dimensions (L×W×H) Weight	$W \times L \times H = 170 \times 150 \times 120 \text{ mm}$ $(6.69" \times 5.91" \times 4.72")^3 \text{ weight}$ = 4 kg (8.82 lb)
Mechanical connection	1/4" – 11/4", others on request
Materials in contact with process gas	Stainless steel AISI 316 L (1.4404), fluoroelastoer, sapphire, PTFE, polymer H.L., nitrite
Filters	PTFE 0.22 μm, PTFE 5 μm
Power supply	24 V 1A

Electronic Connections

Power supply	12-pin M12 male connector
Signal output	12-pin M12 female connector
Electronic output	Active output, maximum 505 Ohm at 24V power supply RS-232, RS-485, USB, Ethernet (with BACCom), Modbus
Storage	0 - +60°C; < 75% RF non-condensing
Maintenance	1-point calibration with fresh air (0.04 vol.% $\rm CO_{2}$, 20.97 vol.% $\rm O_{2}$) once a month (other conditions possible), optional factory calibration once a year
CE FCC ICES	EN61326-1:2006 FCC 15:2009 Subpart 107/109, ICES- 001:2006
Remarks	Don't use in explosive atmosphere, in anoxic atmosphere, in gases with polymers or silicon components or in gases with halogens (F, Cl, Br, etc.), CFC or SOX and H2S

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Sartorius Stedim Biotech GmbH August-Spindler-Str. 11 37079 Goettingen, Germany

Phone +49.551.308.0 Fax +49.551.308.3289 www.sartorius-stedim.com

Copyright by Sartorius Stedim Biotech GmbH. Goettingen, Germany. All rights reserved. No part of this publication may be reprinted or translated in any form or by any means without the prior written permission of Sartorius Stedim Biotech GmbH. The status of the information. specifications and illustrations in this manual is indicated by the date given below. Sartorius Stedim Biotech GmbH reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Status: 12 | 2014, Sartorius Stedim Biotech GmbH, Goettingen, Germany

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Vorort-Reparaturen sind nicht möglich. Informieren Sie den Service der Sartorius Stedim Biotech GmbH.
Lassen Sie die vollständige Kalibrrierung nur vom Service der Sartorius Stedim Biotech GmbH durchführen.